# Instructions for Use

Product: GL-DTRAX® Cervical Cage System

**REF** DX-10-500  
DX-10-600

**Providence Medical Technology, Inc.**  
4234 Hacienda Drive Suite 150  
Pleasanton, CA 94588  
USA  
Tel: +1 (415) 923-9376  
Fax: +1 (415) 923-9377

**Emergo Europe**  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands  
Email: [EmergoVigilance@ul.com](mailto:EmergoVigilance@ul.com)

**QS Engineering AG**  
CH-4106 Therwil  
[info@qs-engineering.ch](mailto:info@qs-engineering.ch)

Before using a product placed on the market by Providence Medical Technology, the operating surgeon should study carefully the following recommendations, warnings and instructions. Providence Medical Technology is not liable for complications that may arise from the use of the device in circumstances outside of the company's control, including, but not limited to, product selection and deviation from the device's indicated uses or surgical technique.

## PACKAGE CONTENTS:

### DX-10-500

Item #	Quantity	Description
1	1	Access Chisel
2	1	Access Chisel Handle
3	1	Decortication Trepine
4	1	Fork Mallet
5	1	Guide Tube
6	1	Rasp Decorticator
7	2	Implant Delivery Tool Subassembly
8	2	GL-DTRAX Cervical Cage-T Implant
9	1	Bone Graft Tamp

### DX-10-600

Item #	Quantity	Description
1	1	Access Chisel
2	1	Access Chisel Handle
3	1	Decortication Trepine
4	1	Fork Mallet
5	1	Guide Tube
6	1	Rasp Decorticator
7	2	Implant Delivery Tool Subassembly
8	2	GL-DTRAX Cervical Cage-SE Implant
9	1	Bone Graft Tamp

## DESCRIPTION:

The GL-DTRAX Cervical Cage System Includes two GL-DTRAX Cervical Cage implants manufactured from implant grade titanium alloy (6Al-4V ELI Titanium). The implants are single-use only. The GL-DTRAX Cervical Cage System also includes a set of sterile, single-use instruments used to access the target spinal level, prepare the joint for fusion, and deliver the GL-DTRAX Cervical Cage implants.

## INTENDED USE:

The GL-DTRAX Cervical Cage System contains posterior fixation implants and delivery instruments that enable distraction of the facet joint, decompression of the nerve root and stabilization of the spine.

## INDICATIONS FOR USE:

The GL-DTRAX Cervical Cage System is indicated for use in skeletally mature patients for posterior cervical treatment at C3-C7 (inclusive) spinal levels for patients with single level radiculopathy due to degenerative disc disease (DDD) as defined by back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies and/or degenerative disease of the facets.

## CONTRAINDICATIONS:

The GL-DTRAX Cervical Cage System should not be used in any of the following instances:

- Infection, local to operative site.
- Patient is allergic to Titanium.
- Local inflammation.
- Fever or leukocytosis.
- Pregnancy.
- Renal Osteodystrophy.
- Cancer of the spine.
- Rheumatoid Arthritis.
- Grossly distorted anatomy due to congenital abnormalities.
- Osteoporosis.

## WARNINGS:

In any surgical procedure, the potential for complications exists. Consideration should be given to the following but not limited to the list below since every patient may respond differently to a cervical fusion procedure.

- Infection or painful, swollen or inflamed implant site.
- Fracture of the implant.
- Loosening, bending, breaking or dislocation of the implant requiring revision surgery.
- Bone resorption.
- Bone over-production.
- Advanced diabetes.
- Obesity.
- Allergic reaction(s) to implant material(s).
- Histological responses possibly involving macrophages and/or fibroblasts.
- Migration of particle wear debris possibly resulting in a bodily response.
- Embolism.
- Subsidence of the device into the vertebral body.
- Bone loss and/or decrease in density due to stress shielding.
- Malalignment of anatomical structures (i.e. loss of normal spinal contours or change in height).
- Physiological reaction to implant devices due to foreign body intolerance including inflammation, local tissue reaction, and possible tumor formation.
- Loss of neurological function by several mechanisms, including direct compression by component parts, stretching of the spinal cord by component parts, vascular spinal cord compromise, or other mechanisms.
- Any medical or surgical condition that would preclude the potential benefit of spinal implant surgery (i.e. elevation of white blood count (WBC) or marked left shift in the WBC differential count).
- Uncooperative patient or patient with neurologic disorders, incapable of following instructions.
- Revision procedures for which an adequate fit of the prosthesis cannot be achieved.
- Metabolic disorders that may impair bone formation.
- Osteomalacia.
- Poor prognosis for good wound healing (e.g., decubitus ulcer, end-stage diabetes, severe protein deficiency and/or malnutrition).
- Adjacent Segment Disease (ASD).

**PRECAUTIONS:**

**PRE-OPERATIVE PRECAUTIONS:**

- The surgeon must evaluate each situation individually based on the patient's clinical presentation in making any decisions regarding implant selection.
- The surgeon must be thoroughly familiar with the implant, instruments, and surgical procedure prior to performing surgery.
- The surgeon should also use medical devices in accordance with their labeled indications and the manufacturer's instructions for use, especially during insertion and removal.
- Patient selection should consider the following factors which may lead to increased risk of failure and can be critical to the eventual success of the procedure: the patient's:
  - a) weight,
  - b) activity level,
  - c) occupation, and
  - d) other patient conditions which may impact on the performance of the GL-DTRAX system.

Implant longevity and stability may be affected by these variables. An overweight patient can produce high loads on the prosthesis, which can lead to failure of the prosthesis. The surgeon must consider the ability and willingness of the patient to follow instructions and to control their weight and activity level.

- Patient should be in the previously described diagnostic categories described under INDICATIONS FOR USE.
- Patient should NOT be in the contraindication groups listed under CONTRAINDICATIONS in this Instructions for Use (IFU).
- Sterilization and handling procedures conforming to accepted standards are mandatory.
- The techniques for implanting this system should be reviewed by the surgeon prior to use of the system.
- The surgeon should inspect the available components prior to surgery to assure that all necessary components are present.
- An adequate inventory of implants should be available at the time of surgery, normally a quantity in excess of what is expected to be used.
- Different metal types should never be used together.
- The surgeon is expected to follow the instructions made available in training manuals and literature relative to implantation of GL-DTRAX Cervical Cage.
- The surgeon is expected to exercise extreme care in the placement of implants, particularly in regard to neural elements.
- Radiographs should be made if there is any question as to the location of the intended or the actual placement of the implants.
- Components should be received and accepted only in packages that have not been damaged or tampered with. Damaged implants and/or instruments should not be used.
- Components must be carefully handled and stored in a manner that prevents scratches, damage and corrosion.

**INTRAOPERATIVE PRECAUTIONS:**

- Extreme caution should be used around the spinal cord and nerve roots. Damage to the nerves will cause loss of neurological functions.
- Breakage, slippage, or misuse of instruments or implant components may cause injury to the patient or operative personnel.
- Imaging system should be used to facilitate surgery.
- Do not mix instruments from different manufacturers. While rare, breakage of instruments may occur especially with extensive use or excessive force.
- It is recommended to place bone graft material in the area to be fused. The graft material should extend from the upper to the lower vertebrae being fused.
- Inspect devices prior to use for damage during shipment or storage or any out-of-box defects that might increase the likelihood of fragmentation during a procedure.
- Proper implant selection must consider design, fixation, patient weight, age, bone quality, size, activity level, preoperative level of health, and also the surgeon's experience and familiarity with the device. Implant longevity and stability may be affected by these variables. Surgeons should inform the patient about these factors.

**POSTOPERATIVE PRECAUTIONS:**

The physician's postoperative directions and warnings to the patient, and the corresponding patient compliance, are extremely important. Detailed instructions on the use and limitations of the device should be given to the patient.

- The patient must be advised of the limitations of the reconstruction and the need for protection of the implant from full weight bearing until adequate fixation and healing have occurred.
- Periodic follow-up, including post-operation X-rays is recommended to monitor the position and state of the implant components, as well as the condition of the bone.
- Periodic post-operative x-rays are recommended for close comparison with early post-op conditions to detect long-term evidence of changes in position, loosening, bending, or cracking of components.
- The patient is expected to follow the detailed instructions, limitations, and warnings from the operating surgeon.
- The patient and the surgeon must understand that the implant is not expected to support the spine if fusion does not occur.
- The risk of bending, loosening or breakage of the implants during postoperative rehabilitation may be increased if the patient is active, if the patient is debilitated, or otherwise unable to use crutches or other such weight supporting devices.
- The patient should not be exposed to mechanical vibrations that may loosen the device. They should also avoid falls or other sudden jolts in spinal position.
- The patient should be warned of possible limitations and instructed to limit and restrict physical activities, especially lifting and twisting motions and any type of sport participation.
- The patient should avoid the consumption of alcohol or nonsteroidal or anti-inflammatory medications or the use of tobacco products during the postoperative healing phase.
- There is a risk of failure of the implant if the fusion of the spine does not occur. It should be recognized that this may occur and is a function of biology. More surgery may be required in such an event. If a nonunion develops or the components loosen, bend, and/or break the device should be removed immediately.
- The surgeon is expected to supply detailed instructions to the patient regarding postoperative activities. The patient should be advised at their inability to bend at the point of spinal fusion and receive training on how to compensate for this loss of motion.
- The potential for multiple complications exist. These are not necessarily due to deficiencies of the implants and may include fracture of the implants due to fatigue, late infection or sensitivity due to fretting-corrosion, prominence of the implants, and displacement of the implants due to failure of the supporting spinal structure.
- Retrieved implants should be properly disposed of and are not to be reused under any circumstance.
- As a precaution, before patients with implants receive any subsequent surgery (such as dental procedures), prophylactic antibiotics may be considered, especially for high risk patients. Any retrieved devices should be treated in such a manner that reuse in another surgical procedure is not possible.
- The patient should immediately consult the surgeon if there is a reoccurrence of radiculopathy or new onset myelopathy.

**MRI WARNING**  
 The GL-DTRAX Cervical Cage System has not been evaluated for safety or compatibility in the MR environment > 3T. MR risks, including heating, migration, and imaging artifacts next to the implant are known, but have not been evaluated for the GL-DTRAX Cervical Cage System.

**INSTRUCTIONS AND TECHNIQUE:**

**CAUTION:** The technique(s) employed in each surgery will always depend upon the medical judgment of the surgeon exercised before and during surgery as to the best mode of treatment for each patient.

The GL-DTRAX Cervical Cage is placed bilaterally into the facet joint using a posterior to anterior trajectory. It is recommended that commercially available bone graft be used post-implant to aid fusion. Bone graft is not supplied as part of the system.

1. Patient Preparation:
  - a. Intubate the patient and put in prone position with the head on a doughnut. Place the patient on a Wilson frame or, alternatively, two rolls. The patient's position, especially the neck, is neutral—not flexed or extended.
  - b. Identify normal spinal landmarks and obtain true anterior-posterior (AP) and lateral fluoroscopic images. These images should line up the facets and endplates in the lateral view and on the AP view and clearly delineate the medial and lateral aspect of the facet joint at the appropriate spinal level.
  - c. Auxiliary tools such as a Steinman pin can be used to mark the skin at the appropriate implantation trajectories or to mark the medial and lateral borders of the facets.
2. Using standard surgical procedure, make an incision through the skin, subcutaneous tissues and fascia. The incision is made approximately a fingerbreadth to two fingerbreadths off the midline and at an appropriate spot to gain trajectory to the facet to be treated.
3. Place the Access Chisel (1) through the ipsilateral skin incision with the chisel blade edge in a vertical orientation in order to get through the musculature. (Figure 1)

Figure 1: Insert Access Chisel

4. **NOTE:** The Access Chisel is designed with a positive stop feature, which is unidirectional. It is intended to engage the inferior articular process of the level above. (Figure 2)

When it is approximated closer to the facet, rotate the tool to orient the chisel blade edge in a horizontal position, parallel to the facet joint gap (Figure 3). Ensure that the black line on the Access Chisel shaft is facing cranially; the black line indicates the side of the Access Chisel with the Positive Stop feature. (Figure 4)

Figure 2: Positive Stop Feature, Access Chisel

Figure 3: Rotation of the Access Chisel


Figure 4: Orientation Indication Line on the Access Chisel

5. Advance the tool until its distal tip is contacting the inferior aspect of the treated level.

**NOTE:** Use AP and lateral fluoroscopy in the approach to the spine to avoid entering the spinal canal and to ensure proper tool placement in the facet joint. (Figure 5)

Figure 5: Proper placement of the Access Chisel

6. Guide the Access Chisel into the facet joint on the ipsilateral side. Manually advance or gently tap the tool into the facet joint using the Fork Mallet tool (4) as a mallet, until the positive stop feature abuts the inferior articular process of the level above. (Figure 6)

Figure 6: Access Chisel advanced into the facet joint

**NOTE:** The entry into the facet joint should involve small measured steps under precise AP and lateral fluoroscopy. The Access Chisel should advance into the facet with controlled hand pressure and should not require significant malleting. The instrument should slide easily into the facet; if it does not, the facet is likely fused, or the tool is not on the correct anatomic position.

7. Remove or detach the Access Chisel Handle by pulling the handle back by the lock sleeve and hold in that position while removing it from the Access Chisel. (Figure 7)

Figure 7: Access Chisel Handle

8. The Decortication Trephine (3) is then advanced over the Access Chisel until the Decortication Trephine (3) contacts the superior and inferior lateral mass. Use AP and Lateral fluoroscopy to confirm correct positioning of the device. The Decortication Trephine is lightly malletted with the Fork Mallet to begin decortication of the facets. This light malleting is then followed with controlled manual rotation of the Decortication Trephine. It is important to maintain downward pressure on the Access Chisel while rotating the Decortication Trephine.
9. Remove the Decortication Trephine.
10. Place the Guide Tube (5) over the Access Chisel, advancing to the facet. Use the technique from Step 3 to aid in a smooth insertion.

- NOTE:** The Guide Tube is designed with a positive stop feature, which is unidirectional. It is intended to engage the inferior articular process of the level above. (Figure 8)

When it is approximated closer to the facet, rotate the tool to orient the Guide Tube forks in a horizontal position, parallel to the facet joint gap. Ensure that the black line on the Guide Tube handle is facing cranially and aligned with the line on the Access Chisel; the black line indicates the side of the Guide Tube with the Positive Stop feature. (Figure 9)

Figure 8: Positive Stop Feature, Guide Tube


Figure 9: Orientation Indication Line on the Guide Tube

- Advance the Guide Tube until its distal tip is contacting the inferior aspect of the treated level.
 

**NOTE:** Use AP and lateral fluoroscopy in the approach to the spine to avoid entering the spinal canal and to ensure proper tool placement in the facet joint.
- Advance the Guide Tube tip into the facet joint to the appropriate AP depth. Gently tap the tool into the facet joint using the Fork Mallet tool as a mallet, until the positive stop feature abuts the inferior articular process of the level above. (Figure 10)

Figure 10: Guide Tube advanced into the facet joint

**NOTE:** The entry into the facet joint should involve small measured steps under fluoroscopy. Malleting is done with small hand movements and hand feel. The instrument should slide easily into the facet; if it does not, the facet is either fused, or the tool is not in the correct anatomic position.

- Verify the correct placement depth in the facet using the positive stop feature and fluoroscopy. Once deployed, the posterior edge of the GL-DTRAX Cervical Cage will be just anterior to the positive stop features. The distal end of the Guide Tube tool indicates the distal tip location of the implant.
 

**NOTE:** The AP placement is of critical importance as the implant function depends on its correct placement in the facet joint.
- Remove the Access Chisel from the Guide Tube.
- Insert the Rasp Decorticator (6) into the Guide Tube completely. The Rasp Decorticator should be removed and reinserted several times to increase the degree of intra-facet decortication. This ensures that the intra-facet surfaces are decorticated by the Rasp Decorticator's surface features.

- Once Guide Tube placement is confirmed and decortication is complete, remove the Rasp Decorticator from the Guide Tube by pulling the tool out, leaving the Guide Tube tool located in the intra-facet space. Confirm via fluoroscopy that the Guide Tube forks are still located in the desired position in the facet.
 

**NOTE:** Removing or separating the Rasp Decorticator and Guide Tube may be tight depending on the patient's anatomy. If difficulty is encountered, use the flat ended side of the Fork Mallet to spread the handles apart (Figure 11), rotating the Fork if necessary.


Figure 11: Separating the Rasp Decorticator and Guide Tube using the Fork Mallet

- Visually verify that the GL-DTRAX Cervical Cage is loaded onto the Implant Delivery Tool correctly; reference the images below. (Figure 12) The images show the Cervical Cage-SE as examples.


Figure 12: Proper loading of GL-DTRAX Cervical Cage on the implant delivery tool

- Prepare the Cage by packing with fusion material, e.g. autogenous bone graft.
- Under AP and lateral fluoroscopic control advance the Cage and delivery instrument into the target joint. Malleting may be required to fully insert implant and distract target joint.
- If malleting is needed maintain downward pressure on the Guide Tube to ensure it remains positioned in the facet joint during malleting.
- The two handles should lock together securely, with the retention snaps actuated.
- Use fluoroscopy to confirm proper placement of the Cage.
- Once proper Cage position is confirmed locate the Cage Release Knob on the handle of the Delivery Instrument and turn the knob counter clockwise to release the Cage from the Delivery Instrument.
- NOTE:** Turn the knob until the knob spins freely. Pull up on the knob to verify that the Cage is completely released.
- Use fluoroscopy to confirm proper placement of the Cage. Remove Implant Delivery Tool from the Guide Tube by squeezing the release buttons on the Implant Delivery Tool handle and pulling the tool out. Deployed Cage will be left in the intra-facet space.
- Prepare the bone graft and insert directly into the Guide Tube. Use the Bone Graft Tamp (9) to advance the bone graft through the Guide Tube and into the facet joint. Make sure to deliver all bone graft material and use AP and Lateral fluoroscopy to confirm proper positioning.
- Remove Bone Graft Tamp from the Guide Tube. Remove the Guide Tube from the facet.
- Repeat steps 3 through 28 on the contralateral facet.
- Close the incision using established surgical procedures.

**STERILITY:**

GL-DTRAX Cervical Cage System is provided **STERILE** and are **SINGLE-USE** devices.

**CAUTION:** The GL-DTRAX Cervical Cage System should never be re-sterilized or reused after contact with body tissues or fluids, but rather should be discarded. Providence does not take any responsibility for the use of implants re-sterilized after contact with body tissues or fluids.

**STORAGE CONDITIONS:**

All implants must be stored in a clean, dry environment and be protected from sunlight and temperature extremes.


**Verify package integrity prior to use. Do not use system components if packaging is damaged. Do not use the GL-DTRAX Cervical Cage System if the expiration date has been exceeded.**

**LIMITED WARRANTY:**

Providence Medical Technology, Inc. products are sold with a limited warranty to the original purchaser against defects in workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.

**SYMBOLS USED ON THIS PRODUCT:**

Reference: ISO 15223-1, *Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied*

Symbol	Symbol Number and Title	Brief Description
	5.1.1 Manufacturer	Indicates the medical device manufacturer.
	5.1.2 Authorized representative in the European Community	Indicates the authorized representative in the European Community.
	5.1.4 Use-By date	Indicates the date after which the medical device is not to be used.
	5.1.5 Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	5.1.6 Catalog Number	Indicates the manufacturer's catalog number so that the medical device can be identified.
	5.2.4 Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.
	5.2.8 Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	5.3.2 Keep away from sunlight	Indicates a medical device that needs protection from light sources.
	5.3.4 Keep Dry	Indicates a medical device that needs to be protected from moisture
	5.4.2 Do not reuse	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	5.4.3 Consult instruction for use	Indicates the need for the user to consult the instructions for use.
	5.4.4 Caution	Indicates the need for the user to consult the instructions for use for important cautionary information.
	CE Mark	CE Mark and Notified Body number

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